


Technical Specifications & Requirements for Software Installation and Use

Cardea 20/20™ ECG System


About Cardea 20/20 ECG System

Cardea 20/20 ECG records and measures a resting ECG from the adult and pediatric (age ≥14) body surface. It provides automatic ECG interpretations which are identified as “Unconfirmed” by the product until they have been over-read and confirmed by a clinician.




ECG Transmitter with lead wires

The patient lead wires terminate in clip/snap connectors that attach to user-selected commercially available ECG electrodes (electrodes are not included in the Cardea 20/20 ECG system).



Software

Cardea 20/20 ECG personal computer software application



Instructions

Operator’s manual and quick-start guide

Product Specifications

Transmitter Frequency Response	0.0 – 150 Hz
Sampling Rate	1000 samples/second
Resolution	± 0.5 microvolt (vs. industry standards of 2.5-5 microvolts)
Real-Time Recording Screen	Displays full ECG trace for live quality control
One-Click Analysis	Stops recording, analyzes ECG, and provides interpretations
ECG Recording Duration	16 seconds
Interactive Editing	Edit global phase picks (P on, P off, Q, S, end of T)
Echocardiographic Data Entry	Can be saved with ECG data

SOFTWARE & HARDWARE REQUIREMENTS

Disk	Minimum 2 GB of free disk space
CPU	Intel® Core™2 Duo CPU @1.50 GHz or greater, 32-bit (x86) or 64-bit (x64) processor, or equivalent
Display	1024×768 or higher resolution
Memory	Minimum 2 GB of system memory
Pointing Device	Windows® compatible mouse
Keyboard	Windows® compatible keyboard
Windows® Operating Systems Supported	Windows® 7, 8 and 10
Ports	1 USB port
Product support	Call (866) 554-3751, Option 2 from 6:00 AM to 5:00 PM (Pacific Time), Monday through Friday. Email: support@cardiacinsightinc.com .

Cardea 20/20™ ECG is a trademark of Cardiac Insight, Inc.

All product names, logos, and brands are property of their respective owners. All company, product and service names used in this document are for identification purposes only. Use of these names, logos, and brands does not imply endorsement. Refer to Cardea 20/20 Operator's Manual and other product documentation for additional information and specifications.

Cardiac Insight reserves the right to make changes in specifications and features shown herein at any time and without notice or obligation.